



eLife®

ORTHOPEDICS

2018

embrace  
life.

**Europe**

- Austria
- Belgium
- Croatia
- Cyprus
- Denmark
- France
- Germany
- Greece
- Hungary
- Italy
- Macedonia
- Netherlands
- Poland
- Russia
- Serbia
- Switzerland
- Spain
- Sweden
- Turkey
- UK

**Asia Pacific**

- Australia
- Japan
- India
- Indonesia
- Korea
- Malaysia
- Singapore
- Taiwan
- Thailand
- New Zealand

**America**

- USA
- Canada
- Mexico
- Uruguay
- Colombia

**Africa**

- Egypt
- Ethiopia

**Middle East**

- Bahrain
- Israel
- Iran
- Jordan
- Kuwait
- Lebanon
- Iraq
- Qatar
- Saudi Arabia
- UAE



**WORLD PARTNERS**



tel: +886-2-2923-1516  
 fax: +886-2-2929-0783  
 info@e-lifebracing.com

**e**life is dedicated to restoring mobility and improving the quality of life for people around the world. Through our commitment to quality, innovation and product safety we deliver life enhancing product solutions that have a positive impact on a broad range of musculoskeletal conditions.

As one of the world's leading providers of orthopedic supports, e-life has built its reputation on innovative product solutions that are founded on a total commitment to quality and safety. We control every stage of the manufacturing process and apply the same exacting standards, from the sourcing of raw materials to research and development, manufacture and testing. Our commitment to quality and safety is paramount and we exceed all international quality and safety standards, including DNV, DIN, CE and FDA

– we are accredited with ISO 9001 and ISO 13485. Our in-depth medical knowledge combined with our advanced product design expertise enables us to deliver functional products for optimal patient care, helping solve mobility problems and deliver a better quality of life.

Orthopedics Supportive and flexible to optimize movement and comfort. We are committed to restoring mobility and improving the quality of life for patients, with a range of braces and supports that have been developed to enhance mobility across all areas of the human body. As a world leader in medicinal orthopedics, we have the expertise and technical knowledge that ensures our products deliver maximum functionality and comfort.



# CONTENT

## SERIES PRODUCTS

Q-Fit Series	P.5
D1 Series	P.11
Bio-Contour Series	P.15



## UNIT

Knee	P.19
Ankle	P.30
Back	P.36
Neck	P.44
Shoulder	P.48
Elbow	P.52
Wrist	P.55
Pediatric	P.60
Foot	P.68



# Q-FIT SERIES



## Q-Fit system

The Q-Fit system is developed by conception of easy and quick donning and doffing. More than that it brings users more benefits. Q-Fit system replaced the traditional lace function so that users can be

convenient to use and save more time and energy during taking on or off the brace. During activities, users are able to enhance effectiveness of performance and excellent compliance

- Glove- like Fit
- Micro-adjustability
- One handed adjustment
- Durable
- Quick on and off
- Powerful closure
- Lightweight

**BOA®**



### How to use the system



- » Push in to engage
- » Turn clock-wisely to tighten
- » Pull up for quick release

## E-WAR021A


### PREMIUM Q-FIT BACK ARMOUR

#### FEATURES

- » The low profile posterior panel, anterior panel and side panels all are connected via a Velcro closure, and can be removed and reattached easily. (E-WAR021A)
- » The Q-Fit system allows you to do quick, one-handed control your comfort for standing or sitting.
- » 3:1 ratio, low-friction Q-Fit System provides powerful, smooth and easily-controlled compression and support.
- » Single dial front closure provides easy loosening and tightening
- » Thin, lightweight, breathable mesh material combines with contouring panels for increased support.



SIZE	CM	BACK PANEL
S	75-85	S
M	85-95	
L	95-105	L
XL	105-115	
XXL	115-125	



## E-WAR021

### Q-FIT BACK ARMOUR




**E-WAR108**  
**Q-FIT WAIST BRACE**

**FEATURES**

- » Front and back anatomically-contoured plastic panels.
- » Pre-shaped panels restrict spinal undesired ROM
- » Easy to adjust intra-abdominal pressure when changing position from sitting.



SIZE	CM
S	56-71
M	69-84
L	81-96
XL	94-109
XXL	107-122




**E-WAR110**  
**Q-FIT LOW PROFILE BACK BRACE**

**FEATURES**

- » Low profile design with 4 stays used to support L1-L5
- » Easy to adjust intra-abdominal pressure when changing position from sitting.
- » Pocket handle for convenient applying



SIZE	CM
S	56-71
M	69-84
L	81-96
XL	94-109
XXL	107-122





**E-WRR052**  
**Q-FIT WRIST SPLINT**

- » Anatomically-contoured aluminum palmar splint
- » Immobilization of the wrist
- » Comfortable, breathable .skin-friendly material

SIZE	CM
S	13-15
M	15-18
L	18-20
XL	20-22

LEFT /RIGHT




**E-WRR058**  
**Q-FIT THUMB AND WRIST SPLINT**

**FEATURES**

- » Anatomically-contoured aluminum palmar splint
- » Immobilization of the wrist and thumb
- » V-notch at fifth MP for palmar spread
- » Breathable open cell foam material

SIZE	CM
S	13-15
M	15-18
L	18-20
XL	20-22

LEFT /RIGHT





**E-ELR006**  
**Q-FIT ELBOW BRACE**

**FEATURES**

- » The SBR pad provides pain relief at the muscle insertion point
- » Extra elastic band fixes brace in position for easy donning and doffing

SIZE	CM
S	23-25.5
M	25.5-28
L	28-30.5
XL	30.5-33







**E-ANR043**  
**Q-FIT ANKLE BRACE**

**FEATURES**

- » Bilateral plastic stays for additional support
- » Extra non-elastic cross-stirrup straps provide maximal stability
- » Made of duration and tough nylon fabric

SIZE	CM
S	15-20
M	20-25
L	25-30
XL	30-35





**E-ANR045**  
**Q-FIT ANKLE STABILIZER**



**FEATURES**

- » Y-Strap stabilizes the ankle and prevent it from rolling inwards by Q-Fit system adjusting tension on the straps.
- » Bilateral anatomically PE plates provide great support
- » Easy rear foot entry allows for quick on-off

SIZE	CM
S	15-20
M	20-25
L	25-30
LEFT /RIGHT	



# D1 Series




**E-KNR005**  
**Q-FIT PATELLAR STRAP**

**FEATURES**

- » Click buckle designed is more convenient for using
- » The rubber tube provides the patellar tendon support to release pain
- » Breathable material for great comfort

SIZE	CM
S	30-34
M	34-38
L	38-42
XL	42-46
XXL	46-50




Rubber tube



## D1 Series

D1 series is an innovative manufacturing process used to produce complicated molded parts from two different materials. D1 series molding optimizes to create powerful combination of hard and soft materials and enhanced tactile and cosmetic product features.

### FEATURES

- » The varus/valgus hinge can adjust lateral and medial shift angle to unload shear force of the knee (E-KN041)
- » Angular hinge allow safety control for knee motion (E-KN041)
- » Medially applied frame corrects for either medical or lateral compartment
- » Lightweight, and ergonomic frame designed for superior knee support
- » Unique buckle design for easy donning and doffing
- » Breathable liner provides optimal texture

SIZE
UNIVERSAL
LEFT/RIGHT

### E-KN040 OA KNEE STANDARD



### E-KN041 OA KNEE SHELTER

### E-WK016

## D1 WALKER

### FEATURES

- » Double-injection design with full shell and fancy style provides comfort for calf area
- » Double pumps provide individual compression to decrease the swelling
- » Shock absorbing sole reduces heel-strike impact force
- » Front panel for more complete protection and support
- » Athletic outsole design provides better shock absorption and improves anti-slip ability
- » 11" and 17" available



### E-WK017

## D1 WALKER WITH AIR CUSHIONED OUTSOLE

- » Air-cushioned outsole design provides excellent shock absorption




SIZE	WOMEN SHOE	MEN SHOE
S	6-8	4-7
M	8-11	7-10
L	11-15	10-13
XL	15+	13+

**E-WR009  
CMC GUIDER**

- » Ergonomic frame provides superior support and excellent comfort of CMC joint
- » Inner metal stay offers extra support and stability of thumb
- » Easy to apply and adjust with simple straps

SIZE	CM
S	16-19
M	19-23
L	23-26

LEFT / RIGHT




**E-AN065  
D1 STIRRUP BRACE**

- » Anatomically-contoured plastic shell
- » Inflatable air bladders allow user to adjust compression pressure to reduce swelling
- » Simple to apply and adjust by two circumferential straps
- » Easy adjustable heel strap

SIZE
UNIVERSAL



**E-AN064  
DYNAMIC ANKLE SHELTER**

- » Adjustable module is designed for customize fit
- » Hinges allow free plantar motion and prevents damaged ankle torsion
- » Bilateral soft pad eliminates joint pressure
- » Moldable foot plate enhances ankle control and comfort

SIZE	WOMEN SHOE	MEN SHOE
S/M	8-11	7-10
L/XL	11-16	10-15



**BIO-CONTOUR  
BRACES**



## BIO-CONTOUR BRACES

Innovative 3D weave technology makes braces provide optimal body fit. Bio-contour series perfectly provides support over the joints, muscles and soft tissue. Bio-contour

braces not only enhance the injury joints and soft tissue recovery ability but also prevent the injury happened during activities.


BIO-CONTOUR BRACES

### E-WA401 BIO-CONTOUR BACK BRACE WITH STRAP

#### FEATURES

- » 3 plastic stays provide extra low back support
- » Bilateral compression bandages provide intra-abdominal pressure
- » Hand pocket design for easy donning and doffing.

SIZE	CM
S	68-78
M	78-88
L	88-98
XL	98-108
XXL	108-118





### E-WA450 BIO-CONTOUR BACK BRACE

#### FEATURES

- » Triangular silicone pad provides great low back support.
- » 3 plastic stays offer extra back stability.
- » Hand pocket design for easy donning and doffing.

SIZE	CM
S	68-78
M	78-88
L	88-98
XL	98-108
XXL	108-118





### E-KN450 BIO-CONTOUR KNEE BRACE

#### FEATURES

- » Bilateral stays support knee ligament
- » Silicone pad maintains patella in position

SIZE	CM
S	33-36
M	36-38
L	38-41
XL	41-46




### E-AN450 BIO-CONTOUR ANKLE BRACE

#### FEATURES

- » Provides great support and compression.
- » Silicone pads provide ankle joint great stability.

SIZE	CM
S	21-23
M	23-25
L	25-27
XL	27-31



BIO-CONTOUR BRACES

## E-SH450 BIO-CONTOUR SHOULDER BRACE

### FEATURES

- » Special silicone pad improves muscle activation
- » Dual adjustable straps for easy using and fit both sides
- » Axillary soft pad design for comfortable wearing



SIZE	CM
S	25-28
M	28-31
L	31-34
XL	34-37



## E-EL450 BIO-CONTOUR ELBOW BRACE

### FEATURES

- » Two silicone pads provide intermittent massage to targeted areas.
- » The 3D knitted material of the elbow brace secures it in place



SIZE	CM
S	15-20
M	20-25
L	25-30
XL	30-35




## E-WR450 BIO-CONTOUR WRIST BRACE

### FEATURES

- » Provide great support and compression.
- » Silicone pad prevents over-stretch on ulnar side



SIZE	CM
S	13-15
M	15-18
L	18-20
XL	20-22




**KNEE**

**KNEE**

## POST OP KNEE BRACE

Post-Operative Knee Brace is designed to provide protected, controlled range-of-motion for patients recovering from knee surgery or those who have knee injuries or instabilities as following indications:




### E-KN099 POST OP KNEE BRACE (NEW)

#### FEATURES

- » Padded lined malleable aluminum cuffs for a perfect fit
- » Pad with anti-slip design prevents brace from sliding and slipping off
- » Streamlined durable aluminum hinge provides maximum protection but with minimum weight
- » "Pull pin design" enables easy release for slide-to-size length adjustments.
- » Simply press and rotate hinge with lock
  - Extension limitation from 0° ~ 105°
  - Flexion limitation from 0° ~ 120°
  - Lock from 0° ~ 90°
  - Adjustable in every 15 degrees



### E-KN099B POST OP KNEE BRACE CUT OFF TYPE



SIZE	LENGTH (CM)	CIRCUM. (CM)
UNIVERSAL		
INFO		
E-KN099	APPROX 48-61	CALF MAX.48 THIGH MAX. 61
E-KN099B	LONG 57 / SHORT 44	

### E-KN097B POST OP KNEE BRACE

#### FEATURES

- » Excellent collateral stabilization via long side-bars prevents valgus or varus deviation.
- » Adjustable ROM by simple setting hinges
- » Telescoping feature for tool-free length adjustment
- » Simply rotate and press-lock settings
  - Extension limitation from 0° ~ 60°
  - Flexion limitation from 0° ~ 120°
  - Immobilization limitation from 0° ~ 60°
  - Adjustable in every 15 degrees

SIZE	LENGTH (CM)	CIRCUM. (CM)
UNIVERSAL	APPROX 48-61	CALF MAX.48 THIGH MAX. 61



### E-KN096 POST OP KNEE BRACE

#### FEATURES

- » Excellent collateral stabilization via long side-bars prevents valgus or varus deviation.
- » Adjustable ROM by simple setting hinges
- » Transcendent overall adjustability
- » "One push button" can be released for length adjustment
- » Simply rotate and press-lock settings
  - Extension limitation from 0° ~ 90°
  - Flexion limitation from 0° ~ 120°
  - Immobilization limitation from 0° ~ 90°
  - Adjustable in every 15 degrees

SIZE	LENGTH (CM)	CIRCUM. (CM)
UNIVERSAL	APPROX 48-61	CALF MAX.48 THIGH MAX. 61

## OA KNEE


In people who have osteoarthritis (OA) knee, the cartilage wears away, leading to an overgrowth of bone underneath. The cartilage becomes rough and then breaks down.

When the cartilage breaks down, it causes pain, swelling, and problems with moving the joint. These symptoms can worsen over time.



Normal knee

Osteoarthritis (OA) knee



3 POINTS OF LEVERAGE

## E-KN036

# UNLIMIT S1 OA KNEE BRACE

### FEATURES

- » Ultra-light weight, single upright with valgus and varus adjustment to unload
- » Laterally applied frame corrects for either medial or lateral compartment
- » 3 points of leverage for offloading.
- » Buckle design eliminates Velcro adjustment when every donning and doffing.
- » One size fits most users and allows ergonomically fitting
- » ROM hinges with plug-in flex/ext stops



SIZE
UNIVERSAL
LEFT/RIGHT

## E-KN042

# DAILY OA KNEE BRACE

(NEW)

### FEATURES

- » Provide offloading of the medial or lateral compartment for mild to moderate osteoarthritis.
- » The single-sided low profile hinge encourages patient compliance
- » "U" Buttress controls the patella in a more natural position.
- » Open anterior thigh section for easy donning and doffing
- » Comfortable, breathable material provides focused compression around the knee joint to address swelling and promote blood flow



SIZE	THIGH(CM)	CALF(CM)
S	30-37	25-30
M	37-44	30-35
L	44-52	35-42
XL	52-60	42-50

LEFT /RIGHT

## E-KN038

# CONQUER OA KNEE BRACE

### FEATURES

- » Dual uprights with valgus and varus adjustment to offload
- » Double upright provides functional support for those patients with ligament instability
- » Support the knee joint and share the loading of affected joint
- » ROM hinges with plug-in flex/ext stops
- » Full shell design provide well stability
- » Could be used at home or for gentle sports

SIZE	THIGH(CM)	CALF(CM)
S	39-47	29-37
M	47-54	37-43
L	54-60	43-50
XL	60-67	50-57


LEFT /RIGHT



## FUNCTIONAL KNEE BRACE

Functional Knee Brace can effectively provide great stability for tears, injuries, or reconstruction of knee ligaments as followings:

- Anterior cruciate ligament (ACL)
- Posterior cruciate ligament (PCL)
- Lateral collateral ligament (LCL)
- Medial collateral ligament (MCL)



### FEATURES

- » Carbon fiber frame with contoured shape and rigid frame for excellent fitting and stability
- » Flexible subshell with adjustable strap for more customized fit
- » Polycentric hinges with plug-in flex/ext stops  
Extension limitation at: 0°, 10°, 20°, 30°, 40°  
Flexion limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°  
Immobilization limitation at: 0°, 10°, 20°, 30°
- » Numbered straps provide desirable counter force and adjust optimal size.
- » Breathable padding liner for superior ventilation and skin contact

## E-KNK035 X1 KNEE BRACE



## E-KN034 FUNCTIONAL KNEE BRACE

### FEATURES

- » 4-leverage design for optimal stabilization
- » Prevent hyperextension via adjustable hinge ROM
- » ROM hinges with plug-in flex/ext stops  
Extension limitation at: 0°, 10°, 20°, 30°, 40°  
Flexion limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°  
Immobilization limitation at: 0°, 10°, 20°, 30°
- » Silicone-grip strap pads for anti-migration
- » Concentric condyle pads correspond to bilateral femoral condyles.

#### 1. Point of leverage

Femur is anchored by a cuff placed at the top of the anterior thigh.

#### 2. Point of leverage


A strap across the back of the lower thigh pushes the femur anteriorly

#### 3. Point of leverage

A strap placed on the anterior tibia applies an active constant load to prevent anterior tibial translation

#### 4. Point of leverage

Tibia is anchored by a cuff placed bottom of the posterior calf



SIZE	THIGH(CM)	CALF(CM)
S	39-47	29-37
M	47-54	37-43
L	54-60	43-50
XL	60-67	50-57

LEFT / RIGHT

## E-KNK03A X1 KNEE BRACE WITH CUP



SIZE	THIGH(CM)	CALF(CM)
S	39-47	29-37
M	47-54	37-43
L	54-60	43-50
XL	60-67	50-57

LEFT / RIGHT




**E-KN002**  
**KNEE BRACE WITH**  
**PATELLAR OPEN**

**FEATURES**

- » 1/8" thick neoprene with nylon lining
- » Increase blood circulation through warm effect
- » Open patellar design
- » Effectively stabilizes the kneecap to eliminate anterior knee pain

SIZE	CM
S	30-34
M	34-38
L	38-42
XL	42-46
XXL	46-50



**E-KN006**  
**PATELLAR TENDON BRACE**

**FEATURES**

- » 1/8" thick neoprene with nylon lining
- » Open patellar design
- » Effectively stabilize the kneecap to eliminate anterior knee pain
- » Wrap around design provides an easy on and off with a custom fit.



**E-KN004**  
**OPEN PATELLAR KNEE**  
**BRACE**

**FEATURES**

- » 1/8" thick neoprene with nylon lining
- » Open patellar design
- » Popliteal cut-out design for comfort fit
- » Wrap around design provides an easy on and off with a custom fit.

SIZE
UNIVERSAL




**E-KN052**  
**HINGED KNEE SUPPORT**

**FEATURES**

- » 3/16" thick neoprene with nylon lining
- » Open padded buttress supports the patella
- » Lightweight hinges hidden in pockets offer superior stability
- » Upper and lower fasteners offer further fit



**E-KN058**  
**WRAP AROUND HINGED**  
**KNEE SUPPORT**

**FEATURES**

- » 3/16" thick neoprene with nylon lining
- » Hinges hidden in pockets offer superior stability
- » Wrap around design provides an easy on and off with a custom fit.
- » Upper and lower fasteners offer further fit



**E-KN005**  
**PATELLAR STRAP**

**FEATURES**


- » Ideal for patellar tendinitis ,jumper knee
- » Reinforce and protect the patellar tendon
- » Simple fitting & time-saving application via hook and loop closure



SIZE
UNIVERSAL



SIZE	CM
S	30-34
M	34-38
L	38-42
XL	42-46
XXL	46-50




Suitable for E-KN006, E-KN052, E-KN058

## E-KN070 STABILIZED KNEE SUPPORT

### FEATURES

- » 1/8" thick neoprene with nylon lining
- » Bilateral spiral steel stays prevent valgus or varus deviation.
- » Open padded buttress supports the patella
- » Upper and lower coronal fasteners offer further fit

SIZE	CM
S	30-34
M	34-38
L	38-42
XL	42-46
XXL	46-50



## E-KN081 KNEE IMMOBILIZER

### FEATURES

- » Three piece design for precise fit
- » Removable medial, lateral and posterior stays
- » Comfortable, washable, skin-friendly terry toweling cloth material
- » Velcro straps with loop closures for simple fitting
- » Additional patellar strapping(E-KN081A)
- » Ideal for post-operative immobilization, conservative treatment of ligament injury




## E-KN081A KNEE IMMOBILIZER WITH PATELLAR STRAP



Additional patellar strapping

SIZE	CM	LENGTH
S	30-34	
M	34-38	12", 14", 16", 18", 20", 22", 24"
L	38-42	
XL	42-46	available
XXL	46-50	




## E-KN201 KNEE BRACE

### FEATURES

- » Made of circular-knitted and four-way stretch elastic material for comfortable fit
- » Comprised of breathable and moisture absorbed material



SIZE	CM
S	28-34
M	34-41
L	41-48
XL	48-55
XXL	55-63



## E-KN250 KNEE SUPPORT WITH SILICONE PATELLAR RING


### FEATURES

- » Made of circular-knitted and four-way stretch elastic material for comfortable fit
- » Padded buttress supports the patella



Pad

SIZE	CM
S	28-34
M	34-41
L	41-48
XL	48-55
XXL	55-63



## E-TH201 THIGH BRACE

### FEATURES

- » Made of circular-knitted, four-way stretch elastic material for comfortable fit
- » Ideal for muscle strain and micro injury of thigh.



SIZE	CM
S	44-49
M	49-54
L	54-59
XL	59-65



## E-SN201 SHIN BRACE

### FEATURES

- » Made of circular-knitted, four-way stretch elastic material for comfortable fit
- » Ideal for muscle strain and micro injury of calf.



SIZE	CM
S	30-34
M	34-38
L	38-42
XL	42-46



## ANKLE



## UNIT - ANKLE

## FEATURES

- » Lightweight and low profile aluminum side struts
- » Rocker bottom design assists gait
- » Decreased heel-strike force on calcaneus
- » Simple to apply and adjust by circumferential straps
- » Dial lock hinges with adjustable ROM limitation (E-WK002)
  - Dorsiflexion limitation at: 0°, 7.5°, 15°, 22.5°, 30°, 37.5°, 45°
  - Plantarflexion limitation at: 0°, 7.5°, 15°, 22.5°, 30°, 37.5°, 45°
  - Immobilization limitation at: 0°, 7.5°, 15°, 22.5°, 30°, 37.5°, 45°
- » Inflatable aircells are available for further support while reducing swelling (E-WK005)
- » Available in 11" & 17" height (E-WK001/ E-WK005)
- » Available in 14" & 17" height (E-WK002)

SIZE	WOMEN SHOE (US)	MEN SHOE (US)
XS	4.5	UP TO 4
S	5-8.5	4.5-7
M	9-11	7.5-10
L	11.5-13.5	10.5-12.5
XL	13.5+	12.5+



E-WK001  
**FIXED WALKER**

WK001-11"

E-WK002  
**ROM WALKER**



Dial lock hinge

E-WK005  
**PNEUMATIC FIXED WALKER**



Pump

**FEATURES**

- » Lightweight and contoured strut design provides better support and fit
- » Front and back panels for superior protection
- » Rocker bottom design assists gait
- » Decreased heel-strike force on calcaneus
- » Simple to apply and adjust by circumferential straps
- » Available in 11" & 17" height
- » Inflatable aircells are available for further support (WK009/ WK009D)
- » Diabetic Walker System includes: (WK009D)
  - Diabetic Walker
  - Insole WK009D
  - Stocking

SIZE	WOMEN SHOE (US)	MEN SHOE (US)
S	6-8	4.5-7
M	8.5-11.5	7.5-10
L	11.5-13.5	10.5-12.5
XL	13.5+	12.5+

E-WK008

**PRO SHIELD WALKER**



E-NS001

**PLANTAR FASCIITIS STRETCH SPLINT**

**FEATURES**

- » Prefabricated splint
- » Adjustable dorsiflexion straps
- » Soft velour-lined foam padding for optimal comfort
- » Minimal house-hold ambulation is possible via pawl surface under the sole of orthosis
- » Indications of use: plantar fasciitis, achilles tendonitis, contractures and heel spurs

SIZE	WOMEN SHOE (US)	MEN SHOE (US)
S	Up to 7.5	Up to 7
M	8-10.5	7.5-10
L	11+	10.5+

E-WK009

**PRO SHIELD PNEUMATIC WALKER**



E-WK009D

**DIABETIC WALKER**



E-NS003

**DORSAL NIGHT SPLINT**

**FEATURES**

- » Speed recovery from plantar fasciitis, achilles tendonitis and overuse injuries
- » Comfortably stretches the calf and soft tissues of the foot
- » Lightweight and breathable soft padding material for comfort
- » Open heel to hold heel in place and prevent migration during the nighttime

SIZE	WOMEN SHOE (US)	MEN SHOE (US)
S / M	6-10	5-9
L / XL	10.5+	9.5-14



**FEATURES**

- » Carbon fiber design is light to provide excellent dynamic motion
- » Easy donning and doffing with Velcro strap (E-AFK001) or With Q-Fit closure system (E-AFK002)
- » Trimmed foot part for fitting most shoes
- » Anatomical shaped posterior strut for excellent fit

**E-AFK001  
X1 CARBON  
AFO**



**E-AFK002  
Q-FIT  
CARBON  
AFO**



SIZE	SHOE SIZE (EU)	HEIGHT(CM)
S	35-37	32
M	37-39	34
L	39-41	36
XL	41-44	38

LEFT /RIGHT

**E-AF002  
AFO LEAF SPRING**

**FEATURES**

- » Made of durable PP material
- » Improve toe clearance in normal gait
- » Designed specifically for right or left

SIZE	SHOE SIZE (EU)	HEIGHT(CM)
S	35-37	30.5
M	37-39	32
L	39-41	34
XL	41-44	39


LEFT /RIGHT

**E-AN043B  
ANKLE BRACE WITH STRAP**

**FEATURES**

- » Made of durable and tough nylon fabric
- » Extra non-elastic cross-stirrup straps provide maximal stability
- » Bilateral plastic stays for additional support
- » Velcro lace closure design for convenient application

SIZE	CM
S	15-20
M	20-25
L	25-30
XL	30-35




**E-AN051  
LACED ANKLE BRACE**

**FEATURES**

- » Made of durable canvas material
- » Anatomically-contoured wide aluminum ankle stays
- » Lace-up design

SIZE	CM
S	15-20
M	20-25
L	25-30
XL	30-35





**E-AN055B  
AIR ANKLE STIRRUP BRACE**

**FEATURES**

- » Inflatable air bladders allow user to adjust compression pressure to reduce swelling
- » Rigid contoured outer shells control inversion or eversion
- » Simple to apply and adjust by two circumferential straps

SIZE
UNIVERSAL





**E-AN060  
AIR / GEL ANKLE BRACE**

**FEATURES**

- » Combined air and gel bladders for cold therapy and even-distributed pressure
- » Simple to apply and adjust by two circumferential straps
- » Easy adjustable heel strap

SIZE
UNIVERSAL





**E-AN201  
ANKLE BRACE**

**FEATURES**

- » Flexible compression provides ankle support and improves proprioception
- » Made of circular-knit, four-way stretch elastic material for comfortable fit
- » Breathable material with moisture absorbed

SIZE	CM
S	15-20
M	20-25
L	25-30
XL	30-35




## ROLLER SPEED SYSTEM

The system is created by conception of saving energy and convenient adjust abdominal pressure. Two sides pulling provides controlled compression and also prevents brace to migration.




### Universal size design

5 sizes in one, users can adjust sizes easily by hook and loop

Suggestion waist circumference:

S	70-85 (cm)
M	85-95 (cm)
L	95-105 (cm)
XL	105-115 (cm)
XXL	115-125 (cm)



### FEATURES

- » Wrap around design with back, abdominal panel and front panel from stomach to clavicle (WA023)
- » Sternal Pad Kit available for maximum flexion control (WA023)
- » Anatomically-shaped rigid panel fits thoracic and lumbar curvature (WA023)


E-WA023

### FRONT & TLSO PANEL



E-WA021

### LSO PANEL



BACK

BACK BRACE

BACK




**E-WA003  
BACK SUPPORT**

**FEATURES**

- » Plastic splints provide extra support
- » Effective vented tension side-pulls provide superior support
- » Made of high quality vented elastic strapping material

SIZE	CM
S	56-71
M	69-84
L	81-96
XL	94-109
XXL	107-122




**E-WA004  
INDUSTRIAL BACK SUPPORT**

**FEATURES**

- » Effective tension side-pulls provide superior support
- » Anti-migrated textured webbing and suspension system
- » Made of highly vented mesh material

SIZE	CM
S	70-81
M	81-96
L	96-119
XL	119-142




**E-WA052  
ABDOMINAL BINDER/  
3 PANELS**

**FEATURES**

- » Provide compression to stabilize abdominal area
- » Made of stretch elastic material
- » Ideal for protection after surgery and support the belly



**E-WA104  
LUMBAR SACRO  
SUPPORT**

**FEATURES**

- » 4 rigid steel rods and 2 plastic rods support back muscles
- » Effective tension side-pulls provide superior support
- » Made of vented mesh material




**E-WA106  
LUMBAR SACRO  
SUPPORT**

**FEATURES**

- » 4 rigid anatomically contoured aluminum rods support back muscles
- » Rhomboid tension side-pulls flatten the lordosis
- » Breathable open cell foam material

SIZE	CM
S	56-71
M	69-84
L	81-96
XL	94-109
XXL	107-122




Suitable for E-WA052, E-WA104, E-WA106




### E-WA110 LUMBAR SACRO SUPPORT

#### FEATURES

- » 4 rigid anatomically contoured aluminum rods support back muscles
- » Effective vented tension side-pulls provide superior support
- » Made of high quality vented elastic strapping material

SIZE	CM
S	56-71
M	69-84
L	81-96
XL	94-109
XXL	107-122




### E-WA113 DORSAL LUMBAR BACK SUPPORT

#### FEATURES

- » Rigid aluminum uprights and pelvic bar limit spine flexion, extension, and lateral bending
- » Bilateral shoulder straps with pads
- » Made of foam padding material for maximal comfort

SIZE	CM
S	56-71
M	69-84
L	81-96
XL	94-109
XXL	107-122




### E-WA120 OSTEOPOROSIS BRACE

#### FEATURES


- » Anatomical aluminum plate fits spinal curve
- » The front open piece makes this brace be easily applied like a vest
- » Elastic bands constrain abnormal movement and posture but remain mobility
- » Ideal for osteoporosis users



Aluminum upright



SIZE	CM
S/M	56-84
L/XL	84-109



### E-HE001B HYPEREXTENSION BRACE

#### FEATURES

- » Rigid frame restricts side-bending movement
- » Moveable pelvic bar for better comfort while sitting
- » Adjustable frame width and height
- » Adjust hyperextension force via "One Push Buckle" (E-HE001)
- » Ideal for stable compression fractures, and osteoporosis




(E-HE001B)



(E-HE001)



SIZE	PELVIC (CM)	LENGTH (CM)
S	60-75	39-46
M	75-90	42-49
L	90-105	45-52
XL	105-115	49-56





### E-WA118 LUMBAR SACRO SUPPORT

#### FEATURES

- » Posterior lumbar panel provides stabilization and comfort
- » Side straps for additional pressure adjustment
- » Soft and breathable material provides comprehensive comfort

SIZE	CM
S	56-71
M	69-84
L	81-96
XL	94-109
XXL	107-122





### E-HP008 HSO HIP SUPPORT

#### FEATURES

- » 2 flexible and 3 rigid, removable aluminum stays help to control hip abduction/adduction
- » Side pockets for convenient change of stays
- » Wrap-around simple design
- » Using of highly permeable fabric and skin-friendly material

SIZE	WAIST (CM)	THIGH (CM)
S	61-81	51-61
M	76-97	56-66
L	91-107	61-71
XL	102-127	66-81

LEFT / RIGHT



### E-HP009 HIP BRACE LINE

#### FEATURES

- » Rigid outer shell provides optimal support
- » 6 blocking bolts limit hip flexion and extension
  - Extension limitation at: 0°, 20°, 40°, 60°, 80°, 100°
  - Flexion limitation at: 0°, 20°, 40°, 60°, 80°, 100°
- » Abduction/ Adduction can be adjusted continuously between -60° to 60°
- » Length adjustments from 42-48 cm
- » Using of highly permeable fabric and skin-friendly material



SIZE
UNIVERSAL

### E-HB002 HERNIA BELT

#### FEATURES

- » Crotch straps are tied under the lower margin of gluteus maximum muscle to reduce excess pressure around buttock
- » Two removable flexible compression cushions
- » Made of velvet fabric material

SIZE	CM
S	63-76
M	76-89
L	89-102
XL	102-115



### E-MB002 DELUXE MATERNITY BELT

#### FEATURES

- » Supports the lumbar and abdominal regions to help relief of low back pain
- » 2 rigid steel rods support back muscles
- » Effective tension side-pulls provide superior support
- » Made from lightweight tension elastic and velvet material

SIZE	GESTATIONAL AGE (WEEKS)
S	3-10
M	11-14
L	16-18
XL	20+



#### FEATURES

- » Provide support to the rib cage
- » Pain relief through compression whenever breathing or thorax movements are painful
- » Made of high quality vented elastic strapping and velvet material
- » Female breast contour design (E-RB002)

### E-RB002 RIB BELT FOR WOMEN



SIZE	E-RB002	E-RB001
S	56-66	68-81
M	66-76	81-96
L	76-86	96-112
XL	86-102	112-127

CHEST CIRCUM. (CM)

### E-RB001 RIB BELT FOR MEN



E-CO007  
**ARMOUR CERVICAL COLLAR (NEW)**

**FEATURES**

- » Two-piece polyethylene design enhances patient comfort and is easily adjusted with buckle design
- » The unique adjustment technology allows 7 levels of chin height adjustment for customized fitting with a simple knob sliding.
- » Removable liner is made of soft, and moisture wicking coolmax for greater comfort ,and proper hygiene.
- » The large tracheal opening provides functional access for the airway management and other essential procedures
- » Safe for use with X-ray, CT, and MRI

**Patent No.**  
**10620692740**

**NECK**




1. Removable foam laminated pads wick away moisture to ensure proper patient hygiene and comfort while preventing skin irritation

2. Buckle design makes applying easier and saves more time. SIZE adjustment is only on the first applying by Velcro tabs

3. The unique adjustment technology allows 7 levels of chin height adjustment for customized fitting with a simple knob sliding



SIZE	NECK CIRCUM.	HEIGHT
UNIVERSAL	30-45 CM	11-17.5 CM




E-CO001/004  
**FOAM CERVICAL COLLAR/  
FIRM DENSITY COLLAR**



**FEATURES**

- » Kinesthetic reminder to limit neck movement
- » Made of foam with a stockinette cover
- » Made of high density foam with a stockinette cover (E-CO004)

SIZE	CM
XS	34-38
S	38-41
M	41-45
L	45-50
XL	50-52




E-CO002A  
**CERVICAL COLLAR WITH  
TRACHEA OPENING**

**FEATURES**

- » Anatomical shape and preformed chin provide patient better comfort
- » An anterior hole for a tracheostomy is available
- » Made of closed cell foam material

SIZE	NECK (CM)	HEIGHT (CM)
S	25-33	
M	33-41	5.7/8.3/
L	41-48	10.8/13.3
XL	48 up	





E-CO003  
**RIGID ADJUSTABLE  
COLLAR**

**FEATURES**

- » Anatomical shape and edge padding for great comfort
- » Ventilation holes increase breathability
- » Made of PE material with rigid support

SIZE	NECK (CM)	HEIGHT (CM)
S	25-30	6
M	30-35	8
L	35-40	9
XL	40-45	10




E-CO006  
**CERVICAL COLLAR**

**FEATURES**

- » One-piece injection molded collar provides proper support
- » Sliding tabs for easy adjustment height
- » With trachea open
- » Made of LDPE material with padding

SIZE
UNIVERSAL



E-CH004  
**CERVICAL THORACIC  
ORTHOSIS**

**FEATURES**

- » Immobilize cervical spine and upper thoracic spine
- » "Pull pin design" enables easy release for height adjustments of chin holder.
- » All straps are adjustable
- » Forehead strap adds stabilization
- » Front and back supports for chin and occiput provide maximum control and suitable comfort

SIZE
UNIVERSAL





## E-AR002 SHOULDER IMMOBILIZER

### FEATURES

- » Padded cotton shoulder strap with velour-lined
- » Envelope style with thumb support
- » Breathable, durable material for arm rest

SIZE	CM	
S	25-28	
M	28-30	
L	30-33	
XL	33-36	



## E-AR003 MESH ARM SLING

### FEATURES

- » Envelope style with thumb support
- » Vented nylon mesh with shoulder strap

SIZE	CM	
S	25-28	
M	28-30	
L	30-33	
XL	33-36	



## E-AR007 SHOULDER SLING

### FEATURES

- » Wide chest swathe for optimal support
- » Cross wearing provides better stabilization
- » Comfortable, washable, skin-friendly terry toweling cloth material

SIZE
UNIVERSAL




## E-AR011 SHOULDER IMMOBILIZER WITH ABDUCTION

### FEATURES

- » Immobilization of shoulder joint at 15° abduction
- » Comfortable foam pillow provides more support
- » An exercise ball facilitates hand function and stimulates circulation
- » Quick release buckle makes it easy don and off
- » Soft velour material

SIZE	CM
S	25-28
M	28-30
L	30-33
XL	33-36



## E-SH006 OVER SHOULDER HUMERAL SPLINT (NEW)

### FEATURES

- » The two pieces of the brace come together to form a cylinder to provide compression.
- » Shell of the deltoid brace is a rigid polyethylene that is lightweight and durable
- » Brace is laminated with closed cell foam lining
- » Deltoid cap extension of the broken shoulder brace prevents it from migrating
- » The strap on proximal splint wraps beneath the opposite arm, limiting the motion of shoulder
- » Ideal for fractured humeral treatment

SIZE	CM
S	25-28
M	28-30
L	30-33
XL	33-36




## E-CL003 CLAVICLE POSTURE SHOULDER BRACE

### FEATURES

- » Help to support an early kyphosis and to improve posture
- » Flexible stays help provide added back support
- » Made of high quality vented elastic strapping material

SIZE	CM
S	61-76
M	76-91
L	91-107
XL	107-122




## E-CL004 CLAVICLE BRACE

### FEATURES

- » Help to alleviate clavicle fractures
- » Excellent aid in correcting poor posture and shoulder slump
- » Made of stretch elastic material for flexible compression

SIZE	CM
S	61-76
M	76-91
L	91-107
XL	107-122



## E-SH001 NEOPRENE SHOULDER SUPPORT

### FEATURES

- » Elevate humeral head for anti-gravity
- » Relief of tensile force at the subacromial space
- » 1/8" thick neoprene with nylon lining
- » Ideal for shoulder dislocation/subluxation


SIZE	CM
S	25-28
M	28-30
L	30-33
XL	33-36



## E-EL094A POST-OP ELBOW BRACE WITH HAND GRIP (NEW)

### FEATURES

- » Provide single brace (without hand grip and strap) for customer's order (E-EL094)



Simply press and rotate hinge with lock

- Extension limitation from 0° - 105°
- Flexion limitation from 0° - 120°
- Lock limitation from 0° - 90°
- Adjustable in every 15 degrees



## E-EL092 ROM ELBOW BRACE

### FEATURES

- » Simply press and rotate settings
  - Extension limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
  - Flexion limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°, 105°, 120°
  - Immobilization limitation at: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°
- » "One push bottom" can be released for length adjustment
- » Breathable laminated soft foam material



SIZE	LENGTH (CM)	CIRCUM. (CM)
UNIVERSAL	34-43	Forearm Max. 32
	LEFT /RIGHT	

Suitable for E-EL092, E-EL094, E-EL094A

**E-EL005**  
**PNEUMATIC ELBOW BRACE**

**FEATURES**

- » The pneumatic pad provides relief at the muscle insertion point
- » Hook and loop closure for pressure adjustment



SIZE	
UNIVERSAL	

**E-EL006**  
**PADDED TENNIS ELBOW BRACE**

**FEATURES**

- » The SBR pad provides compression at the muscle insertion point to relieve pain
- » 1/8" thick neoprene with nylon lining



SIZE	CM	
S	23-25.5	
M	25.5-28	
L	28-30.5	
XL	30.5-33	

**E-EL008**  
**EPI**

**FEATURES**

- » Rigid plastic shell gives excellent support and provides compression
- » Nonwoven fabric pad provides comfort



SIZE	CM	
XS	20-22	
S	22-25	
M	25-28	
L	28-31	
XL	31-34	

**E-EL201**  
**ELBOW BRACE**

**FEATURES**

- » Flexible compression provides elbow support
- » Made of seamless, four-way stretch elastic material for comfortable fit



SIZE	CM	
S	23-26	
M	26-31	
L	31-37	
XL	37-43	



**WRIST**

**WRIST**


**E-WR123**  
**FUNCTIONAL RESTING SPLINT (NEW)**

**FEATURES**

- » Provide functional hand position, immobilization and prevent deformity.
- » Optional finger separators.
- » Prefabricated splint
- » Plastic shell covered with foam padding cushion

SIZE	CM
S	6-7
M	7-8
L	8-9
XL	9-10

LEFT /RIGHT




**E-WR050**  
**WRISTLET WITH ABDUCTED THUMB SLEEVE**

**FEATURES**

- » Immobilization of the thumb
- » Anatomically-contoured aluminum thumb splint
- » Stretch elastic material for excellent compression
- » Fit right and left hand



**E-WR056**  
**WRIST BRACE**

**FEATURES**

- » Immobilization and stabilization of the wrist
- » Anatomically-contoured aluminum palmar splint
- » Stretch elastic material for excellent compression

**E-WR057**  
**WRIST SPLINT**

**FEATURES**

- » Immobilization of the wrist
- » Anatomically-contoured aluminum palmar splint
- » V-notch at fifth MP for palmar spread
- » Breathable open cell foam material



**E-WR002**  
**WRIST BRACE**

**FEATURES**

- » Provide compressive and proprioceptive effect over wrist joint.
- » Wrap around design for precise fit
- » Universal design to fit right and left hand


SIZE
UNIVERSAL

**E-WR007**  
**THUMB BRACE**

**FEATURES**

- » Immobilization of the thumb
- » Anatomically-contoured aluminum thumb splint
- » 1/8" thick neoprene with nylon lining

SIZE	THIGH(CM)
S	13-15
M	15-18
L	18-20
XL	20-22




**E-WR058**  
**THUMB WRIST AND PALM SPLINT**

**FEATURES**

- » Immobilization of the wrist and thumb
- » Anatomically-contoured aluminum palmar splint and thumb stay
- » V-notch at fifth MP for palmar spread
- » Breathable open cell foam material





E-WR060  
**WRIST BRACE WITH THUMB**

**FEATURES**

- » Anatomically-contoured aluminum thumb splint
- » Extra-wide circumferential strap for superior support
- » Breathable material with terry-lined interior

E-WR063  
**WRIST BRACE**

**FEATURES**

- » Anatomically-contoured aluminum palmar splint
- » Extra-wide circumferential strap for superior support
- » Universal design fits right or left hand
- » Breathable material with terry-lined interior



E-WR201  
**WRIST BRACE**

**FEATURES**

- » Stabilize both wrist joint
- » Breathable material with moisture absorbed
- » Made of circular-knit, four-way stretch elastic material for comfortable fit



SIZE	CM
S	13-15
M	15-18
L	18-20
XL	20-22

LEFT /RIGHT

Suitable for :  
WR050/ WR056/ WR057/ WR058/ WR060/ WR063 (For both side)/ WR201 (For both side)



SIZE	CM
S	8.9
M	10
L	12.7

E-FS001  
**BASEBALL FINGER SPLINT**

**FEATURES**

- » Immobilize sprained finger
- » Stabilize the phalange joints in position
- » Made of malleable aluminum with foam padding
- » Packed in 10 pcs/set



SIZE	CM
S	3.8
M	5.0
L	7.6

E-FS004  
**FOUR PRONG FINGER SPLINT**

**FEATURES**

- » Immobilize sprained finger
- » Stabilize the phalange joints in position
- » Made of malleable aluminum with foam padding
- » Packed in 10 pcs/set



SIZE	CM
S	5.4
M	7.0
L	9.5

E-FS002  
**TOAD FINGER SPLINT**

**FEATURES**

- » Immobilize sprained finger
- » Stabilize the phalange joints in position
- » Made of malleable aluminum with foam padding
- » Packed in 10 pcs/set



SIZE	CM
S	6.4
M	7.0
L	8.9

E-FS005  
**FOLD-OVER FINGER SPLINT**

**FEATURES**

- » Immobilize sprained finger
- » Stabilize the phalange joints in position
- » Made of malleable aluminum with foam padding
- » Packed in 10 pcs/set



SIZE	CM
S	6.4
M	10
L	15.2

E-FS003  
**GUTTER FINGER SPLINT**

**FEATURES**

- » Immobilize sprained finger
- » Stabilize the phalange joints in position
- » Made of malleable aluminum with foam padding
- » Packed in 10 pcs/set



SIZE
1 2 3 4 5 5.5 6 7

E-FS006  
**MALLET FINGER SPLINT**

**FEATURES**

- » Stabilize the phalange joints in position
- » Made of high quality rigid plastic.
- » Ventilated design for great comfort.
- » Packed in 5 pcs/set



## E-CO801 SOFT FOAM COLLAR

### FEATURES

- » Kinesthetic reminder to limit neck movement
- » Contour design for optimal fit
- » Made of medium density foam with a stockinette cover

SIZE	NECK CIRCUM.(CM)	HEIGHT(CM)	
PEDIATRIC	21-27	4	
YOUTH	25-31	6	



## E-CO802A CERVICAL COLLAR WITH TRACHEA OPENING

### FEATURES

- » Restrict cervical spine flexion & extension with gentle support
- » Light weight, with trachea opening
- » Kinesthetic reminder to limit neck movement
- » Made of closed cell foam material, prevent shear stress

SIZE	NECK CIRCUM.(CM)	HEIGHT(CM)	
PEDIATRIC	25-33	6	



## E-CO806 ADJUSTABLE COLLAR

### FEATURES

- » Anterior oversized trachea opening for emergency
- » Back openings for ongoing cervical assessment
- » Radiolucent, MRI and CT scan compatible
- » Made of polypropylene material with smooth foam padding

SIZE	NECK CIRCUM.(CM)	
PEDIATRIC	20-46	




### E-HM801 CRANIAL PROTECTION HELMET (NEW)



#### FEATURES

- » Made of material IXPE with excellent impact resistant
- » Article ventilation isolation provides great ventilation between head and hat
- » Y type double loops reduce the migration of the hat during activities
- » Adjustable head circumference is able to fit most
- » Two colors available (blue and pink)

SIZE	HEAD CIRCUM.
UNIVERSAL	43cm~55cm



### E-AR801 ARM SLING WITH PAD FEATURES

- » Soft shoulder strap with cotton pads
- » Envelope style with thumb support
- » Made of comfortable, washable, skin-friendly terry fabric

SIZE	AGE IN YEARS	FOREARM LENGTH (CM)
PEDIATRIC	2-7	13-20
YOUTH	8-12	18-25




### E-CL801 CLAVICLE SUPPORT FEATURES

- » Help to correct and remind poor posture and Immobilize fractured clavicle
- » Using of foam strap padding with stockinet covered

SIZE	AGE IN YEARS	CHEST CIRCUM.
PEDIATRIC	2-7	51-61





### E-EL881 ELBOW IMMOBILIZER (NEW)

#### FEATURES

- » To immobilize or stabilize your child's elbow after an injury, sprain, bone fracture, or surgery.
- » Made from a comfortable, latex free foam material that is soft against the skin.

SIZE	SPLINT HEIGHT (CM)	ELBOW CIRCUM.(CM)
PEDIATRIC	10	Up to 19
YOUTH	15	18-23




### E-WR801 WRIST SPLINT FEATURES

- » Anatomically-contoured aluminum palmar splint
- » Low contour at palmar crease allows full finger function
- » Made of comfortable, washable, skin-friendly terry toweling cloth material

SIZE	AGE IN YEARS	WRIST CIRCUM.(CM)
PEDIATRIC	2-7	11-14
YOUTH	8-12	14-17

LEFT /RIGHT




### E-WR802 WRIST SPLINT WITH ABDUCTED THUMB FEATURES

- » Rigid thumb stay for custom fit
- » Anatomically-contoured aluminum palmar splint
- » Made of comfortable, washable, skin-friendly terry toweling cloth material

SIZE	AGE IN YEARS	WRIST CIRCUM.(CM)
PEDIATRIC	2-7	11-14
YOUTH	8-12	14-17

LEFT /RIGHT




**E-KN851**  
**NEOPRENE HINGED KNEE SUPPORT**

**FEATURES**

- » Open patellar design
- » Hinges hidden in pockets offer superior stability
- » Made of closed cell neoprene



**E-KN870**  
**NEOPRENE STABILIZED KNEE SUPPORT**

**FEATURES**

- » Open patellar design
- » Hinges hidden in pockets offer superior stability
- » Made of closed cell neoprene



SIZE	AGE IN YEARS	KNEE CIRCUM.(CM)
PEDIATRIC	6 - 7	25 - 30
YOUTH	8-12	30 - 33

Suitable for : E-KN851, E-KN870

**E-KN881**  
**KNEE IMMOBILIZER**

**FEATURES**

- » Removable medial, lateral and posterior stays immobilize the leg
- » Available length: 8", 12"
- » Made of comfortable, washable, skin-friendly terry toweling cloth material with pediatric print



**E-KN882**  
**KNEE IMMOBILIZER/ 3 PANELS**

**FEATURES**

- » Three piece design for precise fit
- » Removable medial, lateral and posterior stays
- » Available length: 8", 12"
- » Made of soft velvet material with washable and skin-friendly terry-lined interior



SIZE	LENGTH	KNEE CIRCUM. (CM)
UNIVERSAL	8", 12"	25 - 30

Suitable for : E-KN881, E-KN882

**E-WK802B**  
**WALKING BOOT**

**FEATURES**

- » Contoured strut design provides better support and fit
- » Rocker bottom design assists normal gait
- » Washable foam liner improves fit and comfort
- » Front and back panels for superior protection (E-WK802)

SIZE	CHILD SHOE(US)
S	<5-6
M	7-8
L	9-10
XL	11-12



**E-AN860**  
**AIR / GEL STIRRUP ANKLE BRACE**

**FEATURES**

- » Anatomically-contoured plastic shells control inversion or eversion
- » Combined air and gel bladders for cold therapy and even-distributed pressure
- » Foam padding with tricot liner provides good shock-absorption



SIZE	AGE IN YEARS	ANKLE CIRCUM.
PEDIATRIC	2-7	13-19

**E-CA801**  
**CAST BOOT**

**FEATURES**

- » Rocker bottom design assists normal gait
- » Non-skid sole
- » Durable canvas material with adjustable closures

SIZE	LENGTH(CM)	WIDTH(CM)
PEDIATRIC	18	9
YOUTH	22	10




**E-PH801  
PAVLIK HARNESS**

**FEATURES**

- » Harness is a dynamic flexion abduction orthosis
- » Booties design provides advanced comfort
- » Soft stirrups with contact closure for ease of application and adjustment

SIZE	MONTHS	CHEST(CM)
XS	0-3	30-36
S		36-41
M	3-6	41-46
L	6-9	46-53





**E-PH802  
FREJKA SPLINT (NEW)**

**FEATURES**

- » Ideal for congenital hip luxation and sub-luxation, hip dysplasia.
- » Flexible polyethylene splint with foam-padded edges for comfort covered by nice textile diaper

SIZE	(CM)
S	16
M	18
L	20
XL	22




Tied on both sides of strap


**E-WRW801  
PEDIATRIC WRIST SLEEVE (NEW)**

**FEATURES**

- » Ideal for Slight wrist strain or sprain and protection during activities.
- » Made of Stretch knitted material that is breathable and flexible.



SIZE	WRIST CIRCUM.
S	15-18 cm
L	18-21 cm




**E-ELW801  
PEDIATRIC ELBOW SLEEVE (NEW)**

**FEATURES**

- » Flexible compression provides elbow support
- » Made of Stretch knitted material that is breathable and flexible.



SIZE	ELBOW CIRCUM.
S	14-17 cm
L	17-20 cm




**E-ANW801  
PEDIATRIC ANKLE SLEEVE (NEW)**

**FEATURES**

- » Ideal for general protection of ankle joint during activities
- » Made of Stretch knitted material that is breathable and flexible.



SIZE	ANKLE CIRCUM.
S	14-17 cm
L	17-20 cm




**E-KNW801  
PEDIATRIC KNEE SLEEVE (NEW)**

**FEATURES**

- » Ideal for general protection of knee joint during activities
- » Made of Stretch knitted material that is breathable and flexible.



SIZE	KNEE CIRCUM.
S	20-25 cm
L	25-30 cm





## E-HV001

### BUNION NIGHT SPLINT

#### FEATURES

- » Plastic brace holds your toe in position to stretch out tight tendons and toe muscles.
- » Protect bunion against rubbing
- » Maintain metatarsal-phalangeal joint alignment

SIZE	SHOES (EU)
S	34-36
M	37-40
L	41-46
LEFT /RIGHT	



## E-SI002

### SILICONE HEEL CUP FOR SPURS-CENTRAL

#### FEATURES

- » Cups shape provides heel support
- » Specially designed soft blue dot for additional heel cushion

SIZE	SHOES (EU)
S	34-36
M	37-39
L	40-42
XL	43-45



## E-SI005

### HEEL CUSHION

#### FEATURES

- » Specially designed teardrop plug for additional heel cushion

SIZE	SHOES (EU)
S	34-38
M	38-42
L	42-46



## E-SI007

### "U" SHAPE HEEL PAD

#### FEATURES

- » U-shaped design lifts up the heel and release the pain on the central area
- » Made of medical level silicone rubber

SIZE	SHOES (EU)
S	34-38
M	38-42
L	42-46



**E-SI021**  
**FULL LENGTH SILICONE INSOLE**



- FEATURES**
- » Softer inserts absorb shock and pain on the metatarsal and heel areas.
  - » Raised metatarsal pad helps to support metatarsal arch.

SIZE	XS	S	M	L	XL	XXL	XXXL
SHOES (EU)	34	36	38	40	42	44	46

**E-SI102**  
**HEEL CUP WITH FABRIC**



- FEATURES**
- » Skin-friendly and hygienic fabric helps reduce friction
  - » Specially designed soft blue dot for local pressure and pain relief functions

SIZE	S	M	L	XL
SHOES (EU)	34-36	37-39	40-42	43-45

**E-SI121**  
**FULL LENGTH INSOLE WITH FABRIC**



- FEATURES**
- » Skin-friendly and hygienic fabric helps reduce friction
  - » Specially designed soft blue dot for local pressure and pain relief functions

SIZE	XS	S	M	L	XL	XXL	XXXL
SHOES (EU)	34	36	38	40	42	44	46

**E-FO001**  
**3/4 LENGTH ORTHOTIC ARCH SUPPORTS**



- FEATURES**
- » Pre-molded orthotics
  - » Ultimate medial arch support with mid cushion material PPE
  - » Skin-friendly and non-slip faux suede top surface helps reduce friction

SIZE	SHOES (EU) (Unit: One pair)
	22 24 26 28 30 32 34 36 38 40



**E-TS003**  
**TOE SEPARATOR**

- FEATURES**
- » Made of medical level silicone rubber
  - » Reduce friction between toes

SIZE	S	M	L

**E-TS005**  
**BUNION SHIELD WITH LOOP**

- FEATURES**
- » Protect bunion against rubbing
  - » Maintain metatarsal-phalangeal joint alignment

SIZE	S/M	L/XL



**E-TS009**  
**TUBULAR TOE BANDAGE**

- FEATURES**
- » Absorb shock and shear force to reduce pain
  - » Made of TPR

SIZE	LENGTH
XS	15 CM
S	
M	
L	

**E-TS014**  
**FOOT -BUNION SEPARATOR**

- FEATURES**
- » Absorbs lateral shock and minimizes friction on the bunion area.
  - » It brings constant relief of the affected site.
  - » Recommend S/M size for 5th toe bunion and L/XL for big toe.

SIZE	S/M	L/XL



## **e-life International Co., Ltd.**

**Office:** 7F, No. 1, Baosheng  
Road, Yonghe Dist.,  
New Taipei City 234,  
Taiwan, R.O.C.

**P:** +886 2 2923 1516

**F:** +886 2 2929 0783

**E:** [info@e-lifebracing.com](mailto:info@e-lifebracing.com)



[www.e-lifebracing.com](http://www.e-lifebracing.com)