

2017 SPORTS & MEDICINE

elife



**enhance
life.**

SPORTS Medicine



INTRODUCTION


E-life sport medicine is committed to developing the ideal brace that offers complete protection for athletes and sport lovers. Braces varying in strength are available to offer all-round protection of injured joints, muscles, ligaments, tendons, and other tissues, thereby ensuring excellent stability.

To cater to the needs in sports arenas or outdoor activity environments, e-life sport medicine made of neoprene material produces warmth and moderate pressure on joints and muscles. In addition, cool-fit material has also been developed to provide good ventilation and moisture absorption as well as environmental protection material features. Knitting technology is also adopted to provide moderate pressure and the 3D trim for a more comfortable fit.

Enthusiastically adhering to the R&D spirit and innovative ideas, we strive to create more benefits for sport doers in the field of sports medicine.



www.e-lifebracing.com



Support Level 1	_____	P04
Level 2	_____	P08
Level 3	_____	P13
Knitting Level 1	_____	P20
Level 2	_____	P23
Level 3	_____	P25
Sleeve	_____	P27
Padding & Protection	_____	P30
Cold & Hot Therapy	_____	P33

**enhance
life.**

SUPPORT

Level I
SUPPORT

Stabilization
Dynamic protection
Adjustable compression

- **LEVEL I**
Airprene with air mesh used
Mild support
Excellent ventilation

- **LEVEL II**
XBCOOL™ material used
Moderate support
Great moisture wicking
Keep warm

- **LEVEL III**
Cool-fit eco-friendly material used
Maximal support
High quality
Superior fit

ESSENTIAL KNEE SUPPORT

E-KNN001

- Contoured design for freedom of movement.
- Convenient finger tabs offer easy grip pull when fitting.
- 5 needle safety stitch for durable and stronger seam.

Measure around Knee Joint.

SIZE	S	M	L	XL	XL
CM	30-34	34-38	38-42	42-46	46-50



Breathable airprene provides ventilation and optimal compression.

3D air mesh for breathability and excess heat release.

Level 1
SUPPORT



E-ANN001

Essential Ankle Support

- Contoured design for freedom of movement.

Measure above ankle bone.

SIZE	S	M	L	XL
CM	15-20	20-25	25-30	30-35



E-THN001

Essential Thigh Sleeve

- Convenient finger tabs offer easy grip pull when fitting.

Measure circumference of thigh.

SIZE	S	M	L	XL
CM	44-49	49-54	54-59	59-65



E-SNN001

Essential Calf Sleeve

- Convenient finger tabs offer easy grip pull when fitting

Measure largest circumference of calf.

SIZE	S	M	L	XL
CM	30-34	34-38	38-42	42-46




E-WAN001

Essential Back Support

- Wrap around design for customized compression adjustment.

Measure around waist.

SIZE	S/M	L/XL
CM	68-98	91-127



E-SHN001

Essential Shoulder Wrap

- Adjustable straps for easy using and fit both sides

SIZE	UNIVERSAL
------	-----------



E-ELN001

Essential Elbow Support

- Contoured design for freedom of movement.

Measure largest circumference of forearm.

SIZE	S	M	L	XL
CM	15-20	20-25	25-30	30-35



E-WRN001

Essential Wrist Support

- Wrap around design for customized compression adjustment.

SIZE	UNIVERSAL
------	-----------



XBCOOL™ OPEN KNEE SUPPORT

E-KND051

- Bilateral double spiral stays provide knee great support and stability.
- The wrap around design with adjustable velcro straps makes it easy for people to put on this knee brace and customize the fit.
- Extra straps provide additional compression for knee.

SIZE	S/M	L/XL
CM	30-38	38-46



XBCOOL™ material features breathable, snugly fit and comfortable.



E-KND050

XBCOOL™ Knee Wrap

- Bilateral spiral stays provides knee great support and stability.
- Adjustable alternating strong velcro closures provide even support.
- Open patella design allows maximum compression around the kneecap.

SIZE UNIVERSAL




E-KND005

XBCOOL™ Patellar Strap Classic

- Innovative tubular buttress applies constant compression at the patellar tendon.
- Adjustable hook and loop closure for custom fit.

SIZE UNIVERSAL



E-KND006

XBCOOL™ Patellar Strap

- Innovative tubular buttress applies constant compression at the patellar tendon.
- Adjustable alternating strong velcro closures provide even support.

SIZE UNIVERSAL



SUPPORT
Level 2



E-AND001

XBCOOL™ Ankle Support

- Wrap around design for customized compression adjustment and easy applying.

SIZE UNIVERSAL



E-AND002

XBCOOL™ Ankle Brace

- An adjustable strap allows for a perfect and extra support in the ankle.

Measure above ankle bone.

SIZE	S	M	L	XL
CM	15-20	20-25	25-30	30-35



E-THD001

XBCOOL™ Thigh Support

- Fully adjustable with hook and loop closures for custom fit.

SIZE UNIVERSAL



E-SND001

XBCOOL™ Shin Support

- Fully adjustable with hook and loop closures for custom fit

SIZE UNIVERSAL



E-WAD001

XBCOOL™ Back Support

- Internal panels to enhance lumbar stability.
- Adjustment straps provide customized fit and compression.

Measure around waist.

SIZE	S/M	L/XL
CM	68-98	91-127



E-WAD002

XBCOOL™ Back Brace

- Four plastic stays provide additional lumbar support.
- Dual adjustment straps provide customized fit and compression.

Measure around waist.

SIZE	S/M	L/XL
CM	68-98	91-127



SUPPORT Level 2



E-SHD001

XBCOOL™ Shoulder Wrap

- X-shape strap provides extra compression on shoulder.
- Fit either side.



SIZE UNIVERSAL



E-ELD001

XBCOOL™ Elbow brace

- SBR acupunctured compression pad provides target point pain relief.

SIZE UNIVERSAL



E-ELD002

XBCOOL™ Elbow Strap

- Pad provides the right compression on painful point .
- Extra strap offers additional pressure and prevent brace from slipping.

Measure largest circumference of forearm.

SIZE	S	M	L	XL
CM	15-20	20-25	25-30	30-35



E-ELD003

XBCOOL™ Elbow Wrap

- Anatomically pre-shaped for ideal fit.
- Adjustable straps for personalized comfort and support.

SIZE UNIVERSAL



COOLFIT MAX ANKLE STABILIZER

E-ANC052


- Special cutting design prevents bunching during bending.
- Breathable mesh fabric on tongue and back of ankle for comfort.
- Convenient finger tabs offer easy grip pull when fitting.

Measure above ankle bone.

SIZE	S	M	L	XL	L/R
CM	15-20	20-25	25-30	30-35	



Molded contoured support structure provides enhanced lateral stability.



Cool-fit (eco-friendly material) provides superior ventilation, moisture wicking.

SUPPORT Level 3

SUPPORT
Level 3



E-KNC070

CoolFit Max Knee Enhancer

- Special straps design offers patella stability and provides extra compression on thigh.
- Bilateral hinges provides support without limited mobility.
- Four-way stretch lycra mesh at back of knee for long lasting comfort.
- Open patella design keeps kneecap in place and prevents slippage.

Measure around Knee Joint.

SIZE	S	M	L	XL	XXL
CM	30-34	34-38	38-42	42-46	46-50



E-KNC071

CoolFit Max Knee Brace

- Special straps on upper and lower knee provide extra compression.
- Bilateral polycentric hinges help to support knee instability.
- Open patella design keeps kneecap in place and prevents slippage.
- Four-way stretch lycra mesh at back of knee for long lasting comfort.

Measure around Knee Joint.

SIZE	S	M	L	XL	XXL
CM	30-34	34-38	38-42	42-46	46-50



E-ANC050

CoolFit Max Ankle Brace

- Spring steel stays for serious support.
- Special cutting design prevents bunching during bending.
- Breathable mesh fabric on tongue and back of ankle for comfort.
- Convenient finger tabs offer easy grip pull when fitting

Measure above ankle bone.

SIZE	S	M	L	XL
CM	15-20	20-25	25-30	30-35



E-ANC051

CoolFit Max Ankle Stabilizer

- Figure-8 straps for medial, lateral and arch support.
- Spring steel stays for serious support
- Special cutting design prevents bunching during bending
- Breathable mesh fabric on tongue and back of ankle for comfort
- Convenient finger tabs offer easy grip pull when fitting.

Measure above ankle bone.

SIZE	S	M	L	XL
CM	15-20	20-25	25-30	30-35





E-ANC062

CoolFit Max Ankle Shelter

- Hinge design for normal gait activity.
- Figure-of-eight straps as ankle taping enhances inversion/eversion support.
- Anatomically-contoured supporting shell with hinge joint.
- Specifically for right or left by anatomical design

SIZE	S/M	L / XL	L/R
Women shoe	5-9.5	10-13.5	L/R
Men shoe	5-8.5	9-12+	



E-ANC063

CoolFit Max Active Ankle Brace

- Stirrup design with bilateral semi-rigid inner plates prevents inversion or eversion.
- Simple to apply and adjust by two circumference straps.

SIZE	S/M	L / XL
Women shoe	5-9.5	10-13.5
Men shoe	5-8.5	9-12+



E-WAC050

CoolFit Max Back Brace

- Semi-rigid outer panel provides strong support without ventilation limit.
- Side-pulls enhance adjustable abdominal pressure.

Measure around waist.

SIZE	S/M	L / XL
CM	68-98	91-127



E-WAC051

CoolFit Max Back Enhancer

- Five rigid stays provide great support for low back.
- Special Side-pulls design offers low back and abdomen optimal 360° pressure.

Measure around waist.

SIZE	S/M	L / XL
CM	68-98	91-127



SUPPORT
Level 3



E-WRC050

CoolFit Max Wrist Glove

- Palmar stay provides support and stabilization of wrist.
- Two fastening straps allow you to customize the fit of your brace.

Measure around wrist.

SIZE	S	M	L	XL	L/R
CM	13-15	15-18	18-20	20-22	



E-WRC051

CoolFit Max Wrist Brace

- Palmar stay provides support and stabilization of wrist.
- Hook and loop closures for custom fit and comfort.

Measure around wrist.

SIZE	S	M	L	XL	L/R
CM	13-15	15-18	18-20	20-22	



E-WRC052

CoolFit Max Wrist Splint

- Anatomically-contoured aluminum palmar splint.
- Special cutting design fits right or left hand.
- Extra-wide circumferential strap for superior support.

Measure around wrist.

SIZE	S	M	L	XL
CM	13-15	15-18	18-20	20-22



KNITTING

KNITTING

Level 1



Great Comfort
Contoured Fit
Stylish Design

- **LEVEL I**
Lightweight
Mild support
Even compression
- **LEVEL II**
Keep warm
Extra support
Double compression
- **LEVEL III**
Thin, lightweight
Extra tight for strong support
Gradual compression

TECHKNIT KNEE SUPPORT

E-KNS202

- Contour design provides superior fit.
- Even compression enhances muscle performance.

Measure circumference above knee in 14 cm

SIZE	S	M	L	XL
CM	38-43	41-45	44-49	47-52



Intensive Compression



Made of nylon, spandex with good elasticity, and breathable fabric keeps skin dry, and more comfortable.



E-ANS202

TECHKNIT ANKLE SUPPORT

- Cover Heel design provides heel more support
- Unique knit offers even compression.

Measure above ankle bone.

SIZE	S	M	L	XL
CM	16-21	21-26	26-31	31-37



E-ELS202

TechKnit Elbow Support

- Contour design prevents slippage and bunching
- Made of nylon, spandex with good elasticity and more comfort.

Measure largest circumference of forearm.

SIZE	S	M	L	XL
CM	19-22	22-25	25-28	28-31



E-WRS203

TechKnit Wrist Support

- Made of nylon, spandex with good elasticity and more comfort.

Measure around wrist.

SIZE	S	M	L	XL
CM	13-15	15-17	17-19	17-21



E-HDS201

TechKnit Headband

- Wear it running, at yoga, during your workout to create your own unique look

SIZE UNIVERSAL



FLEXIKNIT KNEE BRACE

E-KNS251

- High quality, super elastic yarn ,4 way stretch guarantees a long life-span
- 3D circular knitting technology provide best fit and comfort.
- Non-slip grip. No itching or irritation
- Metal spiral stays offers support on both sides of knee.

Measure circumference above knee in 14 cm

SIZE	S	M	L	XL
CM	38-43	41-45	44-49	47-52



Pressure gel pad stabilizes patella in right place and provides massage effect.



E-ANS251

FlexiKnit Ankle brace

- Pressure gel pads provide extra support for ankle joint.
- Unique knitting design prevent discomfort and bunching.

Measure above ankle bone.

SIZE	S	M	L	XL
CM	16-21	21-26	26-31	31-37



E-ELS251

FlexiKnit Elbow brace

- Bilateral gel pads provide more compression on both sides of elbow to relief pain.

Measure largest circumference of forearm.

SIZE	S	M	L	XL
CM	19-22	22-25	25-28	28-31



E-WRS205

FlexiKnit wrist Support

- Extra-long wristband provides more support.

Measure around wrist.

SIZE	S	M	L	XL
CM	13-15	15-17	17-19	19-21



E-WRS206

FlexiKnit Palm Support

- Extra strap provides adjustable compression on wrist.
- Thumb strap offers great support around thumb.

Measure around wrist.

SIZE	S	M	L	XL	L/R
CM	13-15	15-17	17-19	19-21	



COMPACTKNIT KNEE SLEEVE

E-KNS601

- Provides true graduated compression to promote blood circulation and oxygen flow. Prevents injuries such as cramping, fatigue, soreness and swelling.
- Targeted compression zones provide knee stability without cushion.
- Special design stabilizes patella in place but unlimited motion during movement.

Measure circumference above knee in 14 cm

SIZE	S	M	L	XL
CM	38-43	41-45	44-49	47-52





E-ELS601

CompactKnit Elbow Sleeve

- Multiple zone pressure provides superior support and makes elbow joint and muscle perform efficiently.
- Ultrathin, breathable sleeve provides user all day comfort.

Measure largest circumference of forearm.

SIZE	S	M	L	XL
CM	19-22	22-25	25-28	28-31



E-SNS601

CompactKnit Calf Sleeve

- Thin and breathable fabric with great and proper compression provides calf excellent stability.
- Targeted compression and support and increases proprioception during movement.

Measure largest circumference of calf.

SIZE	S	M	L	XL
CM	30-34	34-38	38-42	42-46



SLEEVE

- Reduced fatigue**
- Enhance performance**
- Improve recovery**



E-HDA001

Head Band

- Made of Lycra Elastic Stretch, these headbands fit nice and snug to your head.
- Two sides can be wearing.

SIZE UNIVERSAL



E-ELA501

Padded Arm Sleeve

- Extra protection pad protects elbow against impact.
- Elastic band with anti-slip designed on top prevents slipping down.
- Flat seam stitching provides extra comfort and additional strength.

Measure around bicep

SIZE	S / M	L / XL
CM	25-31	31-38



E-ELA002

Compression Arm Sleeve

- Designed to fit your arm and provide ideal compression help you recover quickly and stimulates blood circulation.

Measure around bicep

SIZE	S / M	L / XL
CM	25-31	31-38



E-SNA002

Compression Calf Sleeve

- Made of Lycra Elastic Stretch, it fits nice and snug to your calf.
- Comfortable compression that provides support for calf muscles all day.

Measure around the calf

SIZE	S / M	L / XL
CM	31-36	36-41



E-KNA501

Padded Knee Sleeve

- Extra protection pad protects knee against impact.
- Elastic band with anti-slip designed on top prevents slipping down.
- Flat seam stitching provides extra comfort and additional strength.

Measure around knee joint.

SIZE	S / M	L / XL
CM	30-38	38-46



PADDING & PROTECTION

**Keep Warm
Prevent Injury
Absorb Impact**

Padding & Protection



E-KND560

PowerGuard Knee Brace

- High impact pad for superb impact absorption.
- Elastic band with anti-slip designed for preventing slipping down.
- Neoprene provides thermal/compression therapy.
- Stretch lycra mesh for long lasting comfort.

Measure around knee joint.

SIZE	S	M	L	XL
CM	30-34	34-38	38-42	42-46



E-AND501

PowerGuard Ankle Brace

- Bilateral unique cushioning pads provide ankle joint optimal stability to reduce the chances of a sprain.
- Elastic band with anti-slip designed for preventing slipping down.
- Neoprene provides thermal/compression therapy.

Measure above ankle bone.

SIZE	S	M	L	XL
CM	15-20	20-25	25-30	30-35



E-ELD501

PowerGuard Elbow Brace

- High impact pad for superb impact absorption.
- Elastic band with anti-slip designed for preventing slipping down.
- Neoprene provides thermal/compression therapy.
- Stretch lycra mesh for long lasting comfort.

Measure around forearm.

SIZE	S	M	L	XL
CM	15-20	20-25	25-30	30-35



E-KNK035A

X1 Knee protector

- Constructed from an ultra-lightweight carbon fiber composite with knee cap for excellent protection and great stability

SIZE	S	M	L	XL	L/R
Thigh(cm)	39-47	47-54	54-60	60-67	L/R
Calf (cm)	29-37	37-43	43-50	50-57	



E-KN034S

ACL daily guarder

- Convenient short length, 12 inches
- 4-Points-of-Leverage System provides complete ACL protection

SIZE	S	M	L	XL	L/R
Thigh(cm)	39-47	47-54	54-60	60-67	L/R
Calf (cm)	29-37	37-43	43-50	50-57	



E-MG001

Mouthguard

- Air flow channel allows for maximum breathing.
- Integrated lip guard protects upper and lower lips.

Color	Red	Yellow	Blue
Size	Adult		

Cold & Hot Therapy

- Reduce inflammation
- Increase circulation
- Alleviate pain





E-IP053

Cold & Hot elastic Wrap

- Mesh pouch with hook and loop enclosure keeps reusable cold/hot pack in place.
- Easy to use - place in freezer or heat in microwave.
- Ideal for knee, ankle and neck.

SIZE UNIVERSAL



E-IPM001

Cold Therapy Wrap

- The wrap and the 2 elastic Velcro straps help you keep the ice bag in place easily.
- Special cutting designed is ideal for treatment of multiple area as shoulder, back and knee pain; etc.

SIZE UNIVERSAL



E-IPK001

Knee Ice Wrap

- Knee wrap with long straps provide a custom fit with firm comfortable compression to help relieve swelling and bruising.

SIZE UNIVERSAL



E-IPW001

Back Ice Wrap

- It provides effective icing compression while quickly relieving back pain.

SIZE UNIVERSAL



E-IPS001

Shoulder Ice Wrap

- Designed for shoulder injury specific adjustment and freedom of movement during usage.

SIZE UNIVERSAL



E-IP002

Cold & Hot pad

- The gel ice pack is reusable and it can be frozen, heated or microwaved. Used for head, neck, arm, elbow, wrist, back, knee, shin, ankle etc.

SIZE 18x13cm



E-IP003

Ice Bag

- Superior leak resistant cap with large cap opening
- Ideal for application of cold therapy to help reduce pain and swelling.

SIZE 6" / 9" / 11"





elife International Co. Ltd.

Office: 7F, No. 1, Baosheng
Road, Yonghe Dist.,
New Taipei City 234,
Taiwan, R.O.C.

P: +886 2 2923 1516

F: +886 2 2929 0783

E: info@e-lifebracing.com



www.e-lifebracing.com